FreeStyle™ - Kontrolllösung

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FreeStyle™ - Control Solution

Indication for use:

For use with the FreeStyle, FreeStyle Mini™, or FreeStyle Papillon™ Blood Glucose Monitoring System to verify that the FreeStyle or FreeStyle Mini Meter and FreeStyle Test Strips or the FreeStyle Papillon Meter and FreeStyle Papillon Test Strips are working correctly.

Warnings and Precautions:

Always store the control solution at room temperature (below 30°C or 86°F). Do not freeze. The control solution can be used for up to 3 months after opening the bottle or until the expiration date, whichever comes first. Count forwards 3 months from the date you open a new bottle of control solution. This is your discard date. Write this date on the side of the control solution bottle.

Note: If a test is performed with the control solution HIGH, the monitoring system must be used at a temperature of 15 - 40°C (59 - 104°F). If the operating temperature of a test with the control solution HIGH is below 15°C (59°F), an error message (Er2**) may appear.

For further information about this error message, please consult the Owner's Booklet.

Replace the cap on the control solution bottle and close tightly immediately after use.

Do not dilute the control solution with water or any other liquid.

For self-testing. For use outside of the body only (*in-vitro* diagnostic). Do not take internally. If swallowed, seek medical advice immediately.

Control Solution tests are specified to be accurate only when used between 15 - 40°C (59 - 104°F).

Description:

Glucose control solution allows you to perform quality control tests to verify that your FreeStyle, FreeStyle Mini, or FreeStyle Papillon Blood Glucose Monitoring System is working correctly and your test results are reliable.

For detailed instructions on performing a control solution test, please consult the Owner's Booklet.

The control solution contains a fixed amount of glucose.

When should I perform a Control Solution test?

When you question your results and want to confirm that your monitor and test strips are working properly.

For information on how to obtain Control Solution, call Customer Care.

To perform a control solution test, you need a FreeStyle or FreeStyle Mini Meter, FreeStyle Test Strips and FreeStyle Control Solution or a FreeStyle Papillon Meter, FreeStyle Papillon Test Strips and FreeStyle Control Solution. A control solution test follows exactly the same procedure as a test with a drop of blood. Shake the control solution well before use. Gently dip one sample area only of the test strip into the control solution. Verify that the control solution test result is within the range printed on the test strip vial. If the test result is within

the range printed on the test strip vial, your meter and test strips are working correctly.

How do I know whether the Control Solution test results are accurate?

- If control solution results are out of range, REPEAT THE TEST
- Out of range test results may be caused by:
 - an expired or bad control solution
 - an error in performing the test
 - a malfunction of the meter
 - performing the Control Solution test outside of 15 - 40°C (59 - 104°F)
 - an expired or bad test strip
 - a watered-down control solution
 - the code on test strip vial not matching the code set in the meter

If you have ruled out all these possibilities and the test result is still outside of the specified range, call Customer Service.

Do not use the blood glucose monitoring system to test your blood glucose if the control solution test results are outside of the range printed on the test strip vial.

Contents:

Each 4.0 ml vial of FreeStyle Glucose Control Solution contains: D-Glucose (Low 0.05%; High 0.3%), preservative, FD&C red dye, viscosity-adjusting agent.

Symbole • Description des symboles • Descrizzione dei simboli • Verklaring van de symbolen • Descripción de los símbolos • Description of Symbols:

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Abbott
A Promise for Life

CE 0123
ART05269 Rev. B 08/05

EC REP
Abbott Diabetes Care Ltd.
Range Road
Witney, Oxon
OX29 0YL UK

Abbott Diabetes Care Inc.
Alameda, CA 94502 USA
© 2005 Abbott